



ARK Methotrexate II Control

This ARK Diagnostics, Inc. package insert for the ARK Methotrexate II Control must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of the assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Report any serious incident that has occurred in relation to the device to the manufacturer and the appropriate competent authority as applicable. A Summary of Safety and Performance is available through Eudamed (European database on medical devices), SRN: US-MF-000023925.

Customer Service



48089 Fremont Blvd
Fremont, CA 94538 USA
Tel: 1-877-869-2320
Fax: 1-510-270-6298
customersupport@ark-tdm.com
www.ark-tdm.com



EC REP

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Key to Symbols Used

	Batch code	 YYYY-MM-DD	Use by/Expiration date
	Catalog Number		Manufacturer
	Authorized Representative		CE Mark with notified body number
	In Vitro Diagnostic Medical Device		Temperature limitation
	Consult Instructions for Use		Quality Control
Rx Only	For Prescription Use Only		

© 2025, ARK Diagnostics, Inc

Control Kit 5071-0003-00

Control Kit 5071-0003-01

Control Kit 5071-0003-02

1 Name

ARK *Methotrexate II Control*

2 Intended Use

ARK Methotrexate II Control is intended for use in quality control of the ARK Methotrexate II Assay.

3 Content

ARK Methotrexate II Control is comprised of a synthetic protein matrix with the following concentrations of methotrexate:

REF	Product Description	Quality Control	Quantity/Volume
5071-0003-00 Complete Set	ARK Methotrexate II Control* Methotrexate, buffer, bovine serum albumin, and sodium azide (nominal level)	Expected Range ($\mu\text{mol/L}$)	Dropper vials
5071-0003-01 Calibration Range Controls ¹	LOW (0.070 $\mu\text{mol/L}$)	0.050 – 0.090	1 X 4 mL
	MID (0.400 $\mu\text{mol/L}$)	0.300 – 0.500	1 X 4 mL
	HIGH (0.800 $\mu\text{mol/L}$)	0.600 – 1.000	1 X 4 mL
5071-0003-02 High Range Controls ²	5 $\mu\text{mol/L}$	3.75 – 6.25	1 X 2 mL
	50 $\mu\text{mol/L}$	37.5 – 62.5	1 X 2 mL
	500 $\mu\text{mol/L}$	375 – 625	1 X 2 mL

*To convert results from $\mu\text{mol/L}$ methotrexate to $\mu\text{g/mL}$ methotrexate, divide $\mu\text{mol/L}$ by 2.2005. Methotrexate levels become 0.0318, 0.1818, 0.3636, 2.27, 22.7, and 227 $\mu\text{g/mL}$ respectively.

Each laboratory should establish its own ranges for each new lot of controls.

¹ Calibration range controls may be obtained separately as a set, **REF** 5071-0003-01.

² High range controls may be obtained separately as a set, **REF** 5071-0003-02.

4 Warnings and Precautions

- For *In Vitro* Diagnostic Use.
- Do not mix controls from different lot numbers.
- Use each lot as a set.
- Controls contain $\leq 0.09\%$ sodium azide.

5 Instructions For Use

- For a complete summary and explanation of the ARK Methotrexate II Assay, refer to the package insert for the ARK Methotrexate II Assay, REF 5071-0001-00.
- Controls are ready to use. Mix each level by gentle inversion before dispensing.
- High range controls must be diluted prior to testing.
- Squeeze sufficient volume (~40µL/drop) into individual sample cups for each level. Consult instrument-specific sample volume requirements. Return caps to their original containers and keep tight.
- Store at 2-8°C. Once opened, use within 12 months and prior to the expiration date.

6 Limitations of Procedure

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, calibrators, controls, dilution buffer, storage of product as directed, and good laboratory technique.

7 Trademarks

ARKTM is a trademark of ARK Diagnostics, Inc.
Other brand or product names are trademarks of their respective holders.



ARK Diagnostics, Inc.
Fremont, CA 94538 USA

Revised May 2025
1600-1225-00 Rev 04