



# ARK™ *Levetiracetam II Calibrator*

This ARK Diagnostics, Inc. package insert for the ARK Levetiracetam II Calibrator must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of the assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Report any serious incident that has occurred in relation to the device to the manufacturer and the appropriate competent authority as applicable. A Summary of Safety and Performance is available through Eudamed (European database on medical devices), SRN: US-MF-000023925.


## Customer Service




48089 Fremont Blvd  
 Fremont, CA 94538 USA  
 Tel: 1-877-869-2320  
 Fax: 1-510-270-6298  
 customersupport@ark-tdm.com  
 www.ark-tdm.com



2797



Emergo Europe  
 Westervoortsedijk 60  
 6827 AT Arnhem  
 The Netherlands



MedEnvoy Switzerland  
 Gotthardstrasse 28  
 6302 Zug  
 Switzerland

## Key to Symbols Used

	Batch code	 YYYY-MM-DD	Use by/Expiration date
	Catalog Number		Manufacturer
	Authorized Representative	 2797	CE Mark with notified body number
	In Vitro Diagnostic Medical Device		Temperature limitation
	Consult Instructions for Use		Calibrator
<b>Rx Only</b>	For Prescription Use Only		

## 1 Name

ARK<sup>TM</sup> *Levetiracetam II Calibrator*

## 2 Intended Use

ARK Levetiracetam II Calibrator is intended for use in calibration of the ARK Levetiracetam II Assay.

## 3 Content

ARK Levetiracetam II Calibrator is comprised of a synthetic protein matrix with the following concentrations of levetiracetam:

REF	Product Description	Quantity/Volume	
5070-0002-00	<b>ARK Levetiracetam II Calibrators</b> Levetiracetam, buffer, bovine serum albumin, and preservatives	Dropper vials	
	<b>A</b>	0.0 µg/mL	1 X 4 mL
	<b>B</b>	5.0 µg/mL	1 X 2 mL
	<b>C</b>	12.5 µg/mL	1 X 2 mL
	<b>D</b>	25.0 µg/mL	1 X 2 mL
	<b>E</b>	50.0 µg/mL	1 X 2 mL
	<b>F</b>	100.0 µg/mL	1 X 2 mL

## 4 Standardization

There is no internationally recognized standard for levetiracetam. ARK Levetiracetam II Calibrators are prepared by gravimetric dilution of high purity levetiracetam into a synthetic proteinaceous matrix free of levetiracetam.

## 5 Warnings and Precautions

- For *In Vitro Diagnostic* Use.
- Do not mix calibrators from different lot numbers.
- Use each lot as a set.
- Calibrators contain ≤0.09% sodium azide.

## 6 Instructions For Use

- For a complete summary and explanation of the Levetiracetam II Assay, refer to the package insert for the ARK Levetiracetam II Assay, REF 5070-0001-00 and 5070-0001-01.
- Calibrators are ready to use. Mix each level by gentle inversion before dispensing.
- Squeeze sufficient volume (~40µL/drop) into individual sample cups for each level. Consult instrument-specific sample volume requirements. Return caps to their original containers and keep tight.
- Store vials at 2-8°C. Once opened, use within 12 months and prior to the expiration date.

## 7 Procedure

### Calibration

Perform a full calibration (6- point) procedure; test calibrators in duplicate. Verify the calibration curve with at least two levels of quality controls according to the established laboratory quality assurance plan.

### When to Re-Calibrate

- Whenever a new lot number of reagents is used
- Whenever indicated by quality control results
- Whenever required by standard laboratory protocols

## 8 Limitations of Procedure

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, calibrators, controls, storage of product as directed, and good laboratory technique.

## 9 Trademarks

**ARK**<sup>TM</sup> is a trademark of **ARK** Diagnostics, Inc.

Other brand or product names are trademarks of their respective holders.



**ARK Diagnostics, Inc.**  
**Fremont, CA 94538 USA**

Revised October 2025  
1600-1066-00 Rev 06