

## ARK™ Fentanyl Calibrator

This ARK Diagnostics, Inc. package insert for the ARK Fentanyl Calibrator must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of the assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Report any serious incident that has occurred in relation to the device to the manufacturer and the appropriate competent authority as applicable.

### Customer Service



**ARK Diagnostics, Inc.**

48089 Fremont Blvd  
Fremont, CA 94538 USA

Tel: 1-877-869-2320

Fax: 1-510-270-6298

[customersupport@ark-tdm.com](mailto:customersupport@ark-tdm.com)

[www.ark-tdm.com](http://www.ark-tdm.com)

SRN: US-MF-000023925



**2797**











**EC REP**

Emergo Europe  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

**CH REP**

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

### Key to Symbols Used

	Batch code	 YYYY-MM-DD	Use by/Expiration date
	Catalog Number		Manufacturer
	Authorized Representative		CE Mark with notified body number
	Consult Instructions for Use		Calibrator
	Temperature limitation		In Vitro Diagnostic Medical Device
<b>Rx Only</b>	For Prescription Use Only		

© 2025, ARK Diagnostics, Inc.

Negative Kit  5031-0002-01

Cutoff Kit  5031-0002-02

## 1 Name

### **ARK™ Fentanyl Calibrator**

## 2 Intended Use

The ARK Fentanyl Calibrator is intended for use in calibration of the ARK Fentanyl II Assay.

## 3 Content

The ARK Fentanyl Negative Calibrator is composed of a non-sterile, processed human urine matrix that is negative for fentanyl. The ARK Fentanyl Cutoff Calibrator is composed of non-sterile, processed human urine containing 1.0 ng/mL of fentanyl.

REF	Product Description	Quantity/Volume
5031-0002-01	<b>ARK Fentanyl Calibrator A (Negative)</b> Human urine, stabilizer and sodium azide	Dropper Vials
	<b>Negative</b>	0 ng/mL

REF	Product Description	Quantity/Volume
5031-0002-02	<b>ARK Fentanyl Calibrator B (Cutoff)</b> Fentanyl, human urine, stabilizer and sodium azide	Dropper Vials
	<b>Cutoff</b>	1.0 ng/mL

## 4 Standardization

There is no internationally recognized standard for fentanyl. A certified solution of fentanyl is traceable to HPLC. ARK Fentanyl Calibrators are prepared by volumetric dilution of high purity fentanyl into non-sterile, processed human urine free of fentanyl.

Calibrators are made with non-sterile, processed human urine free of fentanyl. Donors were non-reactive in tests for HIV 1/2, HBsAg, HCV, HIV-1 (NAT), HCV (NAT) and RPR.

## 5 Warnings and Precautions

- For *In Vitro* Diagnostic Use.
- Harmful if swallowed.
- Contains human urine. Handle as potentially infectious.

- This product contains  $\leq 0.09\%$  sodium azide. As a precaution, affected plumbing including instrumentation should be flushed adequately with water to mitigate the potential accumulation of explosive metal azides.

## 6 Instructions For Use

- Refer to the package insert for the ARK Fentanyl II Assay for a complete summary and explanation of the fentanyl assay test system.
- Calibrators are ready to use. Mix each level by gentle inversion before dispensing.
- Squeeze sufficient volume ( $\sim 40\mu\text{L}/\text{drop}$ ) into individual sample cups for each level. Consult instrument-specific sample volume requirements. Return caps to their original containers and keep tight.
- Store at  $2-8^{\circ}\text{C}$ . Once opened, use within 12 months and prior to the expiration date.

## 7 Procedure

### Qualitative Results

Use the 1.0 ng/mL Calibrator B as a Cutoff Calibrator to distinguish negative and positive samples. Run the Low and High Controls as Negative and Positive respectively. Report test results less than the absorbance ( $\Delta A$ ) value for the Cutoff Calibrator as Negative. Report results equal to or greater than the absorbance ( $\Delta A$ ) value for the Cutoff Calibrator as Positive.

### When to Re-Calibrate

- Whenever a new lot number of reagents is used
- Whenever indicated by quality control results
- Whenever required by standard laboratory protocols
- A stored calibration curve was effective up to at least 15 days based on supporting data

## 8 Limitations of Procedure

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, calibrators, controls, storage of product as directed, and good laboratory technique.

## 9 Trademarks

**ARK<sup>TM</sup>** is a trademark of ARK Diagnostics, Inc.

Other brand or product names are trademarks of their respective holders.



ARK Diagnostics, Inc.  
Fremont, CA 94538 USA

Revised March 2025  
1600-0560-00 Rev 06