


ARK™ Xylazine Assay



The ARK Xylazine Assay is a homogeneous enzyme immunoassay intended for the qualitative detection and/or semi-quantitative estimation of xylazine and its metabolites in human urine at a cutoff concentration of 10 ng/mL. The assay is intended for use in laboratories with automated clinical chemistry analyzers.




The ARK Xylazine Assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed positive analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) or Liquid Chromatography/tandem Mass Spectrometry (LC-MS/MS) is the preferred confirmatory method. Confirmatory testing and professional judgment should be exercised with any drug test result, particularly when the preliminary test result is positive.

KEY POINTS

- ✦ Convenient, liquid-stable, read-to-use homogeneous enzyme immunoassay
- ✦ Applicable onboard automated chemistry analyzers
- ✦ 0 – 500 ng/mL semi-quantitative calibration range, 10 ng/mL Cutoff
- ✦ High specificity for Xylazine and 4-Hydroxy Xylazine in human urine
- ✦ Nonhazardous preservatives

Next Generation Assays

QUALITATIVE PRECISION



Qualitative Control Precision vs 10 ng/mL Cutoff Calibrator

SEMI-QUANTITATIVE PRECISION

Human Urine (ng/mL)	Relative % Cutoff	# of Results	Mean (ng/mL)	Semi-quantitative Precision Results
0	-100	160	0	160 Negative
2.5	-75	160	2	160 Negative
5.0	-50	160	4	160 Negative
7.5	-25	160	7	160 Negative
10.0	Cutoff	160	10	118 Negative / 42 Positive
12.5	+25	160	12	1 Negative / 159 Positive
15.0	+50	160	14	160 Positive
17.5	+75	160	17	160 Positive
20.0	+100	160	20	160 Positive

Pooled Urine Samples containing Xylazine were assayed in quadruplicate twice a day for 20 days.

ACCURACY – ANALYTICAL RECOVERY

Concentration Tested (ng/mL)	Mean (ng/mL)	Recovery (%)
6	6	96
9	9	100
12	12	97
20	20	101
60	56	93
150	134	89
400	351	88

METHOD COMPARISON

ARK Xylazine Assay Results	Low Negative (<50% of cutoff concentration by LC-MS/MS (< 5.0 ng/mL))	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration by LC-MS/MS) (5.0 - 9.9 ng/mL)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration by LC-MS/MS) (10.0 - 14.9 ng/mL)	High Positive (Greater than 50% above the cutoff concentration by LC-MS/MS) (≥ 15.0 ng/mL)
Positive	0	1*	2	30
Negative	109	5	0	0

*One sample with an LC-MS/MS value of 7.4 ng/mL tested positive in the ARK Xylazine Assay.

CROSS-REACTIVITY

Metabolites

Compound	Concentration Approximately Equivalent to the Cutoff (ng/mL)	Percent Cross-reactivity (%)
3-Hydroxy xylazine	9	110
4-Hydroxy xylazine	19	53
4-Hydroxy xylazine glucuronide	38	26

α-2 Agonist and α-2 Antagonist Compounds

Compound	Concentration Approximately Equivalent to the Cutoff (ng/mL)	Percent Cross-reactivity (%)
Atipamezole	>100,000	<0.01
Brimonidine	15,190	0.07
Clonidine	2,260	0.44
Detomidine	>100,000	<0.01
Dexmedetomidine	>100,000	<0.01
Etomidate	>100,000	<0.01
Etylone	>100,000	<0.01
Guanfacine	>100,000	<0.01
Medetomidine	>100,000	<0.01
Metamizole	>100,000	<0.01
Romifidine	2,780	0.36
Tizanidine	9,350	0.11
Tolazoline	>100,000	<0.01
Yohimbine	>100,000	<0.01

Opiates/Structurally Similar Compounds

Compound	Concentration Tested (ng/mL)	Result (POS/NEG)	Percent Cross-reactivity (%)
6-Acetyl morphine	100,000	NEG	<0.01
Buprenorphine	100,000	NEG	<0.01
Codeine	150,000	NEG	<0.01
Dextromethorphan	100,000	NEG	<0.01
EDDP	100,000	NEG	<0.01
EMDP	100,000	NEG	<0.01
Ethyl morphine	100,000	NEG	<0.01
Fentanyl	100,000	NEG	<0.01
Heroin	100,000	NEG	<0.01
Levorphanol	100,000	NEG	<0.01
Meperidine	100,000	NEG	<0.01
Methadone	100,000	NEG	<0.01
Morphine	100,000	NEG	<0.01
Nalbuphine	100,000	NEG	<0.01
Naloxegol	100,000	NEG	<0.01
Naloxone	100,000	NEG	<0.01
Naltrexone	100,000	NEG	<0.01
Norbuprenorphine	100,000	NEG	<0.01
Norfentanyl	100,000	NEG	<0.01
Norcodeine	100,000	NEG	<0.01
Normorphine	100,000	NEG	<0.01
Noroxycodone	100,000	NEG	<0.01
Nortilidine	100,000	NEG	<0.01
Oxycodone	100,000	NEG	<0.01
Oxymorphone	100,000	NEG	<0.01
Pentazocine	100,000	NEG	<0.01
Tapentadol	100,000	NEG	<0.01
Thebaine	100,000	NEG	<0.01
Tilidine	100,000	NEG	<0.01
Tramadol	100,000	NEG	<0.01

SAFETY AND STABILITY

Reagent on-board stability

Up to at least 60 days

Shelf Life of Reagents, Calibrators, and Controls

18 months from date of manufacturing

Safety

Nonhazardous preservatives

Contains sodium azide $\leq 0.09\%$

Results shown are typical and may vary among laboratory analyzers.

ORDERING INFORMATION

For Criminal Justice and Forensic Use Only within the United States

ARK™ Xylazine Assay (R1 x 28 mL / R2 x 14 mL)	5088-0001-00
ARK™ Xylazine Assay (R1 x 115 mL / R2 x 58 mL)	5088-0001-01
ARK™ Xylazine Calibrator (5 x 10 mL)	5088-0002-00
ARK™ Xylazine Calibrator (2 x 10 mL; Negative)	5088-0002-01
ARK™ Xylazine Calibrator (2 x 10mL; 10 ng/mL Cutoff)	5088-0002-02
ARK™ Xylazine Controls (2 x 10mL; LOW 5 ng/mL) (2 x 10mL; HIGH 15 ng/mL)	5088-0003-00

ARK Diagnostics, Inc.
48089 Fremont Boulevard
Fremont, CA 94538
Tel: 510-270-6270 Fax: 510-270-6298

For Customer Support:
Call toll free: 877-869-2320
customersupport@ark-tdm.com
www.ark-tdm.com

